

GOAT ANTI-LPPR3 ANTIBODY

SKU: EB10566

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.**Unit Size** 100 µg**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.**Synonym /****Alias Names** plasticity-related gene 2 protein| PAP-2-like protein 2| PRG2| PRG-2| LPR3| -|LPPR3**Accession ID** NP_079164.1**Blocking Peptide** EBP10566**Immunogen** Peptide with sequence C-TQRGHDSVYQQN, from the internal region of the protein sequence according to NP_079164.1.**Peptide Sequence** C-TQRGHDSVYQQN**Purification Method** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.**Shipping Instructions** Refrigerated**Predicted Species** Human, Rat**Human Gene ID** 79948**Rat Gene ID** 314614**Product Grade** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png**ELISA****Detection Limit** Antibody detection limit dilution 1:32000.**Application Type** Pep-ELISA

DOCUMENTS

- [Data Sheet](#)