



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-LPPR3 ANTIBODY

SKU: EB10566

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias plasticity-related gene 2 protein| PAP-2-like protein 2| PRG2| PRG-2| LPR3| -|LPPR3

Names

Accession NP 079164.1 ID

Blocking

EBP10566 **Peptide**

Immunogen | NP_079164.1. Peptide with sequence C-TQRGHDSVYQQN, from the internal region of the protein sequence according to

Peptide C-TQRGHDSVYQQN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Rat

Human

79948 **Gene ID**

Rat Gene ID 314614

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA

Type