

GOAT ANTI-LSM2 (AA33-46) ANTIBODY

SKU: EB11571

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	YBL026W U6 snRNA-associated Sm-like protein LSM2 snRNP core Sm-like protein Sm-x5 snRNP small nuclear ribonuclear protein D homolog protein G7b OTTHUMP00000029159 LSM2 homolog, U6 small nuclear RNA associated (S. cerevisiae) LSM2 G7b C6orf28
Accession ID	NP_067000.1
Blocking Peptide	EBP11571
Immunogen	Peptide with sequence C-DQYLNKLTDISVT, from the internal region of the protein sequence according to NP_067000.1.
Peptide Sequence	C-DQYLNKLTDISVT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Reactive Species	Human, Mouse
Human Gene ID	57819
Mouse Gene ID	27756
Rat Gene ID	684148
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 10kDa band observed in lysates of cell line K562 and NIH3T3 (calculated MW of 10kDa according to NP_067000.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

