

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-LSM2 (AA33-46) ANTIBODY

SKU: EB11571



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

YBL026W|U6 snRNA-associated Sm-like protein LSm2|snRNP core Sm-like protein Sm-x5|snRNP|small nuclear Synonym / ribonuclear protein D homolog|protein G7b|OTTHUMP0000029159|LSM2 homolog, U6 small nuclear RNA

Alias associated (S. cerevisiae)|LSM2|G7b|C6orf28 Names

Accession

NP_067000.1 ID

Blocking

EBP11571 **Peptide**

Peptide with sequence C-DQYLNIKLTDISVT, from the internal region of the protein sequence according to **Immunogen**

NP 067000.1.

Peptide

C-DQYLNIKLTDISVT Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Dog, Pig, Cow

Reactive

Species

Human, Mouse

Human

57819 Gene ID

Mouse

27756

Gene ID

Rat Gene ID 684148

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Approx 10kDa band observed in lysates of cell line K562 and NIH3T3 (calculated MW of 10kDa according to

Grade ELISA

Detection Antibody detection limit dilution 1:2000.

Western

Limit

NP_067000.1). Recommended concentration: 1-3µg/ml. Blot

Application

Pep-ELISA, WB

Type



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa