



Telephone: (650) 697-3600



GOAT ANTI-LSM3 ANTIBODY

SKU: EB11572

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

YLR438C|USS2|U6 snRNA-associated Sm-like protein LSm3|SMX4|OTTHUMP00000160275|LSM3 homolog, U6 small

nuclear RNA associated (S. cerevisiae)|LSM3

Names

Accession NP 055278.1

Blocking

Peptide

EBP11572

Peptide with sequence C-SLDERIYVKMRND, from the internal region of the protein sequence according to Immunogen

NP 055278.1.

Peptide

C-SLDERIYVKMRND Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Dog, Pig, Cow

Human

27258 Gene ID

Mouse

Gene ID

67678

Rat Gene ID 297455

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection

Antibody detection limit dilution 1:128000.

Limit

Application Pep-ELISA

Type