

GOAT ANTI-LSM3 ANTIBODY

SKU: EB11572



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	YLR438C USS2 U6 snRNA-associated Sm-like protein LSm3 SMX4 OTTHUMP00000160275 LSM3 homolog, U6 small nuclear RNA associated (S. cerevisiae) LSM3
Accession ID	NP_055278.1
Blocking Peptide	EBP11572
Immunogen	Peptide with sequence C-SLDERIYVKMRND, from the internal region of the protein sequence according to NP_055278.1.
Peptide Sequence	C-SLDERIYVKMRND
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Human Gene ID	27258
Mouse Gene ID	67678
Rat Gene ID	297455
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)