

GOAT ANTI-LTBP2 ANTIBODY

SKU: EB10852



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	MSTP031 LTBP3 C14orf141 latent transforming growth factor beta binding protein 2 LTBP2
Accession ID	NP_000419.1
Blocking Peptide	EBP10852
Immunogen	Peptide with sequence C-NSTSSTEDLPDHD, from the internal region of the protein sequence according to NP_000419.1.
Peptide Sequence	C-NSTSSTEDLPDHD
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Human Gene ID	4053
Mouse Gene ID	16997
Rat Gene ID	59106
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Preliminary experiments in Mouse and Rat Eye lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)