

GOAT ANTI-LUBRICIN / PRG4 ANTIBODY

SKU: EB08496



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	proteoglycan 4, (megakaryocyte stimulating factor, articular superficial zone protein, camptodactyly, arthropathy, coxa vara, pericarditis syndrome) megakaryocyte stimulating factor lubricin camptodactyly, arthropathy, coxa vara, pericarditis syndrome bG174L6.2 (MSF: megakaryocyte stimulating factor) articular superficial zone protein Jacobs camptodactyly-arthropathy-pericarditis syndrome bG174L6.2 SZP MSF JCAP HAPO FLJ32635 CACP proteoglycan 4 PRG4
Accession ID	NP_005798.3; NP_001121180.2; NP_001290161.1
Blocking Peptide	EBP08496
Immunogen	Peptide with sequence DYESFCAEVHNP, from the internal region of the protein sequence according to NP_005798.3; NP_001121180.2; NP_001290161.1.
Peptide Sequence	DYESFCAEVHNP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	10216
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)