

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-LUBRICIN / PRG4 ANTIBODY**

**SKU:** EB08496

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

proteoglycan 4, (megakaryocyte stimulating factor, articular superficial zone protein, camptodactyly, arthropathy,

Synonym / coxa vara, pericarditis syndrome)|megakaryocyte stimulating factor|lubricin|camptodactyly, arthropathy, coxa Alias vara, pericarditis syndrome|bG174L6.2 (MSF: megakaryocyte stimulating factor )|articular superficial zone

**Names** protein|Jacobs camptodactyly-arthropathy-pericarditis

syndrome|bG174L6.2|SZP|MSF|JCAP|HAPO|FLJ32635|CACP|proteoglycan 4|PRG4

Accession

NP\_005798.3; NP\_001121180.2; NP\_001290161.1

**Blocking** 

EBP08496 **Peptide** 

Peptide with sequence DYESFCAEVHNP, from the internal region of the protein sequence according to **Immunogen** 

NP\_005798.3; NP\_001121180.2; NP\_001290161.1.

**Peptide** 

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human **Species** 

Human

10216

**Gene ID Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:128000.

Limit

Application

Pep-ELISA **Type**