



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-LYN ANTIBODY**

**SKU:** EB08450



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / LYN|v-yes-1 Yamaguchi sarcoma viral related oncogene homolog|FL|26625|JTK8|Yamaguchi sarcoma viral (v-Alias

yes-1) related oncogene homolog Names

**Accession** 

ID

NP\_002341.1

**Blocking** 

EBP08450 **Peptide** 

Peptide with sequence C-KTQPVRNTERT, from the internal region of the protein sequence according to **Immunogen** 

NP 002341.1.

**Peptide** 

Sequence

C-KTQPVRNTERT

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Cow

**Species** 

Reactive Human **Species** 

Human

4067

Gene ID **Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:16000.

Limit

Western Approx 60kDa band observed in lysates of cell line U937 (calculated MW of 58.6kDa according to

**Blot** NP\_002341.1). Recommended concentration: 0.1-0.3µg/ml.

**Application** 

Type

Pep-ELISA, WB

## **GALLERY IMAGES**



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa