

GOAT ANTI-LYN ANTIBODY

SKU: EB08450



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	LYN v-yes-1 Yamaguchi sarcoma viral related oncogene homolog FLJ26625 JTK8 Yamaguchi sarcoma viral (v-yes-1) related oncogene homolog
Accession ID	NP_002341.1
Blocking Peptide	EBP08450
Immunogen	Peptide with sequence C-KTQPVRNTERT , from the internal region of the protein sequence according to NP_002341.1.
Peptide Sequence	C-KTQPVRNTERT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Reactive Species	Human
Human Gene ID	4067
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 60kDa band observed in lysates of cell line U937 (calculated MW of 58.6kDa according to NP_002341.1). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

