

# GOAT ANTI-LYPLAL1 ANTIBODY

**SKU:** EB08406



250kDa  
150kDa  
100kDa  
75kDa

50kDa

37kDa

25kDa

20kDa

## SPECIFICATIONS

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.

**Synonym /**

**Alias** KIAA1238|lysophospholipase-like 1|LYPLAL1

**Names**

**Accession ID** NP\_620149.1

**Blocking Peptide** EBP08406

**Immunogen** Peptide with sequence DRFKITNDCEHLE, from the internal region of the protein sequence according to NP\_620149.1.

**Peptide Sequence** DRFKITNDCEHLE

**Purification Method** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

**Shipping Instructions** Refrigerated

**Predicted Species** Human

**Reactive Species** Human

**Human Gene ID** 127018

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

**ELISA**

**Detection Limit** Antibody detection limit dilution 1:32000.

**Western Blot** Approx. 27kDa band observed in Human Adipose, Heart and Kidney lysates (calculated MW of 26.3kDa according to NP\_620149.1). Recommended concentration: 1.5-3.0µg/ml. Primary incubation was 1 hour.

**Application Type** Pep-ELISA, WB

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

