



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-LYPLAL1 ANTIBODY**

**SKU:** EB08406



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / **Alias** 

KIAA1238|Iysophospholipase-like 1|LYPLAL1

Names

ID

**Accession** 

NP\_620149.1

**Blocking** 

EBP08406 **Peptide** 

Peptide with sequence DRFKITNDCPEHLE, from the internal region of the protein sequence according to

NP 620149.1.

**Peptide** 

**Immunogen** 

DRFKITNDCPEHLE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human **Species** 

Reactive

Human **Species** 

Human

127018 **Gene ID** 

**Product** 

Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:32000.

Limit

Western Approx. 27kDa band observed in Human Adipose, Heart and Kidney lysates (calculated MW of 26.3kDa **Blot** according to NP\_620149.1). Recommended concentration: 1.5-3.0µg/ml. Primary incubation was 1 hour.

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**Application** 

Pep-ELISA, WB **Type** 

## **GALLERY IMAGES**



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa