



Telephone: (650) 697-3600

GOAT ANTI-MAMMAGLOBIN B / SCGB2A1 ANTIBODY

SKU: EB06808

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

mammaglobin B|mammaglobin 2|lipophilin C|UGB3|MGC71973|MGB2|LPHC|HGNC:7051|secretoglobin, family 2A,

Alias member 1|SCGB2A1

Names

Accession NP 002398.1 ID

Blocking

EBP06808 **Peptide**

Peptide with sequence CFLNQSHRTLKN, from the internal region of the protein sequence according to Immunogen

NP 002398.1.

Peptide

CFLNQSHRTLKN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human Gene ID

4246

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:2000.

Limit

Application Pep-ELISA **Type**