

GOAT ANTI-MAMMAGLOBIN B / SCGB2A1 ANTIBODY

SKU: EB06808



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	mammaglobin B mammaglobin 2 lipophilin C UGB3 MGC71973 MGB2 LPHC HGNC:7051 secretoglobin, family 2A, member 1 SCGB2A1
Accession ID	NP_002398.1
Blocking Peptide	EBP06808
Immunogen	Peptide with sequence CFLNQSHRTLKN, from the internal region of the protein sequence according to NP_002398.1.
Peptide Sequence	CFLNQSHRTLKN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	4246
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)