

GOAT ANTI-MAPCS (AA97-108) ANTIBODY

SKU: EB11199



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Apcs serum amyloid P component Sap APCS amyloid P-component serum protein
Accession ID	NP_035448.2
Blocking Peptide	EBP11199
Immunogen	Peptide with sequence C-QSKVTVRGMEYY, from the internal region of the protein sequence according to NP_035448.2.
Product Comments	This antibody is expected not to cross-react with CRP and mucosal pentraxin.
Peptide Sequence	C-QSKVTVRGMEYY
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Mouse
Mouse Gene ID	20219
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Preliminary experiments in Mouse and Rat Kidney and Liver lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)