


GOAT ANTI-MAPRE3 (AA151-164) ANTIBODY

SKU: EB09828



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names MAPRE3|RP3|microtubule-associated protein, RP/EB family, member 3|EBF3-S|EBF3|EB3|APC binding protein

Accession ID NP_036458.2

Blocking Peptide EBP09828

Immunogen Peptide with sequence C-KLIGTAVPQRTSPT, from the internal region of the protein sequence according to NP_036458.2.

Peptide Sequence C-KLIGTAVPQRTSPT

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Dog, Cow

Reactive Species Human, Mouse, Rat

Human Gene ID 22924

Mouse Gene ID 100732

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:16000.

Western Blot Approx 33kDa band observed in Human (Cerebellum), Mouse and Rat Brain lysates (calculated MW of 32.0kDa according to Human NP_036458.2, Mouse NP_579928.1 and Rat NP_001007657.1). Recommended concentration: 0.1-0.3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

