

GOAT ANTI-MARCH2 / MARCH-II ANTIBODY

SKU: EB08156



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	membrane-associated RING-CH protein II RNFL172 MARCH-II HSPC240 membrane-associated ring finger (C3HC4) 2 MARCH2
Accession ID	NP_001005415.1; NP_001005416.1; NP_057580.3
Blocking Peptide	EBP08156
Immunogen	Peptide with sequence C-KIREADSPEGPQH, from the internal region of the protein sequence according to NP_001005415.1; NP_001005416.1; NP_057580.3.
Peptide Sequence	C-KIREADSPEGPQH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	51257
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Preliminary experiments in Human Heart, Skeletal Muscle, Lung and Spleen lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)