

GOAT ANTI-MARCH5 / MARCH-V ANTIBODY

SKU: EB08154



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	ring finger protein 153 mitochondrial ubiquitin ligase membrane-associated RING-CH protein
Names	V RNF153 MITOL MARCH-V FLJ20445 RP11-543N17.3 membrane-associated ring finger (C3HC4) 5 MARCH5
Accession ID	NP_060294.1
Blocking Peptide	EBP08154
Immunogen	Peptide with sequence C-RLWRKYSNKLQ, from the internal region of the protein sequence according to NP_060294.1.
Peptide Sequence	C-RLWRKYSNKLQ
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	54708
Mouse Gene ID	69104
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx. 30kDa band observed in Human Heart lysates (calculated MW of 31.2kDa according to NP_060294.1). Recommended concentration: 0.1-0.3µg/ml. An additional band of 200kDa was consistently observed, however this band was not blocked by the immunizing peptide and it is therefore a non-specific signal. We call for caution when used for other assays than Western blot.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

