

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-MARK2 / PAR-1 ANTIBODY

SKU: EB08426





Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

serine/threonine kinase|protein-serine/threonine kinase|ELKL motif kinase 1|ELKL motif **Alias** kinase|MGC99619|EMK1|MAP/microtubule affinity-regulating kinase 2|PAR-1|MARK2 Names

Accession

NP_004945.4; NP_001156768.1; NP_001156769.1 ID

Blocking

EBP08426 **Peptide**

Peptide with sequence C-QNGKDSTAPQR, from the internal region of the protein sequence according to **Immunogen**

NP_004945.4; NP_001156768.1; NP_001156769.1.

Product This antibody is expected to recognise isoform c (NP_004945.4), isoform e (NP_001156768.1) and isoform f

Comments (NP_001156769.1).

Peptide

C-QNGKDSTAPQR Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Cow **Species**

Reactive

Human **Species**

Human

2011

Gene ID

Mouse 13728 Gene ID

Rat Gene ID 60328

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:16000.

Limit

Western Approx 70kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 77.6kDa according to

NP_004945.3). Recommended concentration: 0.3-1µg/ml. **Blot**

Application

Pep-ELISA, WB **Type**



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa