

GOAT ANTI-MARVELD3 (ISOFORM 2) ANTIBODY

SKU: EB11948



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MARVEL domain-containing protein 3 MARVEL (membrane-associating) domain containing 3 MRVLDC3 MARVD3 MARVEL domain containing 3 MARVELD3
Accession ID	NP_443090.4
Blocking Peptide	EBP11948
Immunogen	Peptide with sequence C-DGADGEKAQQLD, from the internal region of the protein sequence according to NP_443090.4.
Product Comments	This antibody is expected to recognize reported isoform 2 (NP_443090.4) only.
Peptide Sequence	C-DGADGEKAQQLD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Pig
Human Gene ID	91862
Mouse Gene ID	73608
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)