

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-MAST205 ANTIBODY

SKU: EB06221

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias microtubule associated testis specific serine/threonine protein kinase|KIAA0807|MTSSK|MAST205

Names

Accession NP 055927 ID

Blocking

EBP06221 **Peptide**

Immunogen Peptide with sequence C-ARAGYGSLSVHKQ, from the C Terminus of the protein sequence according to NP 055927.

Peptide C-ARAGYGSLSVHKQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human

Species

Human

23139 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA

Detection Limit

Antibody detection limit dilution 1:8000.

Western **Blot**

Western Blot: Preliminary experiments gave no signal but low background in Human Testis lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?