

GOAT ANTI-MAST205 ANTIBODY

SKU: EB06221



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	microtubule associated testis specific serine/threonine protein kinase KIAA0807 MTSSK MAST205
Accession ID	NP_055927
Blocking Peptide	EBP06221
Immunogen	Peptide with sequence C-ARAGYGSLSVHKQ, from the C Terminus of the protein sequence according to NP_055927.
Peptide Sequence	C-ARAGYGSLSVHKQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	23139
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Western Blot: Preliminary experiments gave no signal but low background in Human Testis lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?