

GOAT ANTI-GOAT / MBOAT4 ANTIBODY

SKU: EB09665



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	MBOAT4 membrane bound O-acyltransferase domain containing 4 FKSG89 GOAT OACT4 O-acyltransferase (membrane bound) domain containing 4 ghrelin O-acyltransferase
Accession ID	NP_001094386.1
Blocking Peptide	EBP09665
Immunogen	Peptide with sequence C-HRISVFSRKWNQST, from the internal region of the protein sequence according to NP_001094386.1.
Peptide Sequence	C-HRISVFSRKWNQST
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	619373
Mouse Gene ID	234155
Rat Gene ID	306515
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

