

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GOAT / MBOAT4 ANTIBODY

SKU: EB09665

SPECIFICATIONS

Unit Size 100 μg

Synonym / Alias

MBOAT4|membrane bound O-acyltransferase domain containing 4|FKSG89|GOAT|OACT4|O-acyltransferase

(membrane bound) domain containing 4|ghrelin O-acyltransferase **Names**

Accession

ID

NP_001094386.1

Blocking

Peptide

EBP09665

Peptide with sequence C-HRISVFSRKWNQST, from the internal region of the protein sequence according to

Peptide

Sequence

C-HRISVFSRKWNQST

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat

Antibody detection limit dilution 1:2000.

Human

Gene ID

619373

Mouse Gene ID

234155

Rat Gene ID 306515

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA Detection Limit

Type

Application Pep-ELISA