


# GOAT ANTI-MCPT4 / MYONASE (MOUSE) ANTIBODY

**SKU:** EB09192



## SPECIFICATIONS

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.

**Synonym /**

**Alias** myonase|Mcp-4|MMCP-4B|MMCP-4A|MMCP-4|mast cell protease 4|Mcpt4

**Names**

**Accession ID** NP\_034909.2

**Blocking Peptide** EBP09192

**Immunogen** Peptide with sequence C-KLQKKAKETPSVN, from the internal region of the protein sequence according to NP\_034909.2.

**Peptide Sequence** C-KLQKKAKETPSVN

**Purification Method** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

**Shipping Instructions** Refrigerated

**Predicted Species** Mouse

**Reactive Species** Mouse

|                         |                                                                                                                                                                                                                                   |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Mouse Gene ID</b>    | 17227                                                                                                                                                                                                                             |
| <b>Product Grade</b>    | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |
| <b>IHC Results</b>      | Paraffin embedded Mouse Intestine. Recommended concentration: 6-8µg/ml.                                                                                                                                                           |
| <b>ELISA</b>            |                                                                                                                                                                                                                                   |
| <b>Detection Limit</b>  | Antibody detection limit dilution 1:2000.                                                                                                                                                                                         |
| <b>Application Type</b> | Pep-ELISA, IHC                                                                                                                                                                                                                    |

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

