

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-MDC1 / NFBD1 ANTIBODY

SKU: EB06280

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias nuclear factor with BRCT domians protein 1|KIAA0170 gene product|KIAA0170|MDC1|NFBD1

Names

Peptide

Accession NP 055456 ID

Blocking EBP06280

Immunogen Peptide with sequence C-EAFVLSPLEMSST, from the C Terminus of the protein sequence according to NP 055456.

Peptide C-EAFVLSPLEMSST Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Refrigerated

Instructions

Predicted Human **Species**

Human 9656

Gene ID **Product**

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade **ELISA**

Detection Antibody detection limit dilution 1:4000. Limit

Western Blot: Preliminary experiments gave no signal but low background in Human Brain and H460 cell lysates at Western up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with Blot

other antibodies/lysates?

Application Pep-ELISA **Type**