

GOAT ANTI-MDC1 / NFBD1 ANTIBODY

SKU: EB06280



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	nuclear factor with BRCT domains protein 1 KIAA0170 gene product KIAA0170 MDC1 NFBD1
Accession ID	NP_055456
Blocking Peptide	EBP06280
Immunogen	Peptide with sequence C-EAFVLSPLEMSST, from the C Terminus of the protein sequence according to NP_055456.
Peptide Sequence	C-EAFVLSPLEMSST
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	9656
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Western Blot: Preliminary experiments gave no signal but low background in Human Brain and H460 cell lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA