

GOAT ANTI-MED12 (AA1981-1994) ANTIBODY

SKU: EB10987



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	human opposite pai activator-recruited cofactor 240 kDa component OTTHUMP00000023493 OTTHUMP00000023490 OPA-containing protein CAG repeat protein 45 TRAP230 OPA1 OKS KIAA0192 HOPA FGS1 CAGH45 ARC240 mediator complex subunit 12 MED12
Accession ID	NP_005111.2
Blocking Peptide	EBP10987
Immunogen	Peptide with sequence C-NPTLVDPTRHLQQR, from the internal region (near C Terminus) of the protein sequence according to NP_005111.2.
Peptide Sequence	C-NPTLVDPTRHLQQR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Human Gene ID	9968
Mouse Gene ID	59024
Rat Gene ID	679693
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)