


GOAT ANTI-MEIS1 ANTIBODY

SKU: EB09039



250kDa
150kDa
100kDa
75kDa

50kDa
37kDa

25kDa
20kDa

15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names leukemogenic homolog protein|WUGSC:H_NH0444B04.1|Meis1, myeloid ecotropic viral integration site 1 homolog|MGC43380|Meis homeobox 1|MEIS1

Accession ID NP_002389.1

Blocking Peptide EBP09039

Immunogen Peptide with sequence C-SEDITRSANLTDQ, from the internal region of the protein sequence according to NP_002389.1.

Peptide Sequence C-SEDITRSANLTDQ

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Cow

Reactive Species Human

Human Gene ID 4211

Mouse Gene ID 17268

Rat Gene ID 686117

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:32000.

Western Blot Approx. 40kDa band observed in Human Brain (Cerebral Cortex) lysates (calculated MW of 43.0kDa according to NP_002389.1). Recommended concentration: 0.1-0.3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

