

GOAT ANTI-MITF ANTIBODY

SKU: EB05411



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	WS2A[microphthalmia-associated transcription factor]MITF
Accession ID	only full length in Swissprot see QC
Blocking Peptide	EBP05411
Immunogen	Peptide with sequence C-QSESGIVPDFEVG, from the C Terminus of the protein sequence according to only full length in Swissprot see QC.
Peptide Sequence	C-QSESGIVPDFEVG
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Mouse
Human Gene ID	4286
Mouse Gene ID	17342
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western blot: Approx 60kDa band seen in NSO lysate. Recommended for use at 0.5-2µg/ml.

