

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-MLANA ANTIBODY**

**SKU:** EB05002

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias MART-1|MART1|melan-A|MLANA

Names

Accession NP 005502

Immunogen Peptide with sequence C-KLSAEQSPPPYSP, from the C Terminus of the protein sequence according to NP\_005502.

**Peptide** C-KLSAEQSPPPYSP

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Instructions Refrigerated

**Predicted** Human

**Species** 

Human

2315 **Gene ID** 

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade

Western

Western blot: No signal obtained yet but low background observed in Jurkat and HeLa cell lysates at upto 3µg/ml. Blot

Application Pep-ELISA

**Type**