

GOAT ANTI-MLANA ANTIBODY

SKU: EB05002



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MART-1 MART1 melan-A MLANA
Accession ID	NP_005502
Immunogen	Peptide with sequence C-KLSAEQSPPPYSP, from the C Terminus of the protein sequence according to NP_005502.
Peptide Sequence	C-KLSAEQSPPPYSP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	2315
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Jurkat and HeLa cell lysates at upto 3µg/ml.
Application Type	Pep-ELISA