

GOAT ANTI-MOUSE IGG ANTIBODY

SKU: EB2ND-52-3001



SPECIFICATIONS

Formulation	Phosphate buffer saline, PH 7.2, with 0.02% sodium azide as preservative10 mM Tris, PH 8.2, with 150 mM NaCl, 10 mg/ml BSA and 0.1% sodium azide. Quantity: 0.5mg in 0.5mg/ml.
Unit Size	500 µg
Storage Instructions	Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles.
Usage Summary	Suitable for immunoblot (western or dotblot), ELISA, and other peroxidase-based immunoassays. ELISA: 1:16000 WB: 1:3000 to 1:10000. Please optimize dilutions for your particular assay conditions. Starting test dilution: 0.5µg/ml for ELISA and WB.
Immunogen	Purified whole mouse IgG (unconjugated).
Purification Method	Ammonium sulphate precipitation of serum followed by antigen affinity purification.
Shipping Instructions	Refrigerated
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png