

GOAT ANTI-MPP2 ANTIBODY

SKU: EB05252

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MAGUK p55 subfamily member 2 discs large, homolog 2 DKFZp761D0712 DLG2 membrane protein, palmitoylated 2 (MAGUK p55 subfamily member 2) MPP2
Accession ID	NP_005365
Blocking Peptide	EBP05252
Immunogen	Peptide with sequence PVAATNSETAMQQ-C, from the N Terminus of the protein sequence according to NP_005365.
Peptide Sequence	PVAATNSETAMQQ-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	4355
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: Two consistent bands observed in HeLa at approx 70kDa and approx 50kDa with 12 hours incubation of primary.