


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-MRP8 / ABCC11 ANTIBODY

SKU: EB08705



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / multi-resistance protein 8|ATP-binding cassette, sub-family C, member 11|ATP-binding cassette protein

C11|WW|MRP8|EWWD|ATP-binding cassette, sub-family C (CFTR/MRP), member 11|ATP-binding cassette transporter Alias

Names MRP8|ABCC11

Immunofluorescence: Strong expression of the protein seen in the membranes and cytoplasm Usage of HeLa cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric **Summary**

analysis of HeLa cells. Recommended concentration: 10ug/ml.

Accession

NP_115972.2; NP_660187.1; NP_001357425.1 ID

Blocking Peptide

EBP08705

Peptide with sequence C-KVVEFDRPEVLRK, from the internal region (near C Terminus) of the protein sequence **Immunogen**

according to NP_115972.2; NP_660187.1; NP_001357425.1.





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Peptide

C-KVVEFDRPEVLRK

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method

the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted Species

Human

Reactive

Human **Species**

Human

85320 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png Grade

IHC Results Paraffin embedded Human Breast. Recommended concentration: 3.75µg/ml.

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Blot

Approx 125kDa band observed in Human Testis lysates and in preliminary testing of Breast lysate (calculated MW of Western

150kDa according to NP_660187.1). This molecular weight is routinely observed by other sources. Recommended

concentration: $0.5-1.5\mu g/ml$. Primary incubation 1 hour at room temperature.

Application


Pep-ELISA, WB, IHC, IF, FC Type

GALLERY IMAGES



Telephone: (650) 697-3600





250kDa 150kDa 100kDa 75kDa 50kDa

Alexa Fluor 488

37kDa

25kDa

20kDa

15kDa

