

GOAT ANTI-MS4A3 ANTIBODY

SKU: EB05202



SPECIFICATIONS

Unit Size	100 µg
Synonym /	hematopoietic cell 4 transmembrane protein IgE receptor beta subunit IgE receptor beta chain CD20 antigen
Alias	homolog CD20 antigen-like CD20L HTM4 membrane-spanning 4-domains, subfamily A, member 3 (hematopoietic
Names	cell-specific) MS4A3
Accession ID	AAA62319
Immunogen	Peptide with sequence CNSREEISSPPNSV, from the C Terminus of the protein sequence according to AAA62319.
Peptide Sequence	CNSREEISSPPNSV
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	932
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: Band observed in Jurkat at approx 30kDa with 12 hours incubation of primary.