

GOAT ANTI-MS4A7 (CFFM4) ANTIBODY

SKU: EB05139



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	high affinity immunoglobulin epsilon receptor beta subunit CD20/Fc-epsilon-RI-beta family member 4 CD20L4 MS4A8 CFFM4 membrane-spanning 4-domains, subfamily A, member 7 MS4A7
Accession ID	NM_021201
Immunogen	Peptide with sequence CIQQVKSSRSWI, from the C Terminus of the protein sequence according to NM_021201.
Peptide Sequence	CIQQVKSSRSWI
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	58475
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 28-30kDa band seen in Human Testis lysate. Recommended for use at 2-4µg/ml.