

GOAT ANTI-MSH4 (C TERMINUS) ANTIBODY

SKU: EB06348



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	mutS homolog 4 RP11-262K24.2 mutS (E. coli) homolog 4 mutS homolog 4 (E. coli) MSH4
Accession ID	NP_002431.2
Blocking Peptide	EBP06348
Immunogen	Peptide with sequence C-FPRTEQVPEKTEE, from the C Terminus of the protein sequence according to NP_002431.2.
Peptide Sequence	C-FPRTEQVPEKTEE
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	4438
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments gave no signal but low background in human testis lysates up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

