

Telephone: (650) 697-3600

GOAT ANTI-MSH4 (INTERNAL) ANTIBODY

SKU: EB08334

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

e erest

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias mutS homolog 4|RP11-262K24.2|mutS homolog 4 (E. coli)|MSH4

Names

Accession NP 002431.2 ID

Blocking

EBP08334 **Peptide**

Peptide with sequence C-ETNSSTFMKEMKE, from the internal region of the protein sequence according to Immunogen

NP 002431.2.

Peptide

C-ETNSSTFMKEMKE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse

Human

4438 Gene ID

Mouse

55993

Gene ID

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

Application Pep-ELISA **Type**

> Goat Anti-MSH4 (internal) Antibody https://everestbiotech.com/products/goat-anti-msh4-internal-antibody/