

GOAT ANTI-MSH4 (INTERNAL) ANTIBODY

SKU: EB08334



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	mutS homolog 4 RP11-262K24.2 mutS homolog 4 (E. coli) MSH4
Accession ID	NP_002431.2
Blocking Peptide	EBP08334
Immunogen	Peptide with sequence C-ETNSSTFMKEMKE, from the internal region of the protein sequence according to NP_002431.2.
Peptide Sequence	C-ETNSSTFMKEMKE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	4438
Mouse Gene ID	55993
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)