

GOAT ANTI-MTA2 / MTA1L1 ANTIBODY

SKU: EB05164



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	metastasis-associated 1-like 1 metastasis associated gene family, member 2 DKFZp686F2281 metastasis associated gene 1-like 1 metastasis-associated protein 2 MTA1-L1 MTA2 PID metastasis-associated 1-like 1 MTA1L1
Accession ID	NP_004730.2
Immunogen	Peptide with sequence C-LKVKPTLIAVRP, from the C Terminus of the protein sequence according to NP_004730.2.
Peptide Sequence	C-LKVKPTLIAVRP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	9219
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Western blot: Approx 75-85 kDa band seen in Jurkat and Hela Lysate. Recommended for use at 2-4 µg/ml.
Application Type	Pep-ELISA, WB