

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-MTA3 ANTIBODY

SKU: EB06540

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias MTA3-pending|metastasis associated family, member 3|KIAA1266|MTA3

Names

Accession NP_065795 ID

Blocking

EBP06540

Peptide

Peptide with sequence YAAIRAECKMLLNS, from the C Terminus of the protein sequence according to **Immunogen**

NP 065795.

Peptide

YAAIRAECKMLLNS Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Human, Mouse, Rat

Species Reactive

Human **Species**

Human

57504

Gene ID **Product**

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western Blot: Approx 55kDa band observed in Human Breast lysates (calculated MW of 58.8kDa according to Western

Blot NP 065795). Recommended concentration: 1-3µg/ml.