



GOAT ANTI-MTA3 ANTIBODY

SKU: EB06540



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MTA3-pending metastasis associated family, member 3 KIAA1266 MTA3
Accession ID	NP_065795
Blocking Peptide	EBP06540
Immunogen	Peptide with sequence YAAIRAECKMLLNS, from the C Terminus of the protein sequence according to NP_065795.
Peptide Sequence	YAAIRAECKMLLNS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	57504
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western Blot: Approx 55kDa band observed in Human Breast lysates (calculated MW of 58.8kDa according to NP_065795). Recommended concentration: 1-3µg/ml.