

GOAT ANTI-MTMR1 (ISOFORM 1) ANTIBODY

SKU: EB05238



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	myotubularin-related protein 1 myotubularin related protein 1 MTMR1
Accession ID	NP_003819.1
Blocking Peptide	EBP05238
Immunogen	Peptide with sequence C-SSPSHSATSVHTSV, from the C Terminus of the protein sequence according to NP_003819.1.
Product Comments	Please note this will only recognise Human isoform 1 (NP_003819.1), not isoform 2 (NP_789746).
Peptide Sequence	C-SSPSHSATSVHTSV
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	8776
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)