

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-MTMR1 (ISOFORM 1) ANTIBODY

SKU: EB05238

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias myotubularin-related protein 1|myotubularin related protein 1|MTMR1

Names

Accession NP 003819.1 ID

Blocking

EBP05238 **Peptide**

Peptide with sequence C-SSPSHSATSVHTSV, from the C Terminus of the protein sequence according to Immunogen

NP 003819.1.

Product Comments

Please note this will only recognise Human isoform 1 (NP_003819.1), not isoform 2 (NP_789746).

Peptide

Sequence

C-SSPSHSATSVHTSV

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human Species

Human

8776 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA **Type**