

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-MTMR12 ANTIBODY

SKU: EB05237

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

phosphatidylinositol-3-phosphate associated protein|myotubularin related protein 12|3-phosphatase adapter Synonym /

subunit|KIAA1682 protein|KIAA1682|FLJ20476|PIP3AP|3-PAP|phosphatidylinositol-3 phosphate 3-phosphatase **Alias**

Names adaptor subunit|3PAP|MTMR12 protein

Accession

NP 061934.2 ID

Peptide with sequence C-LAKREDEFVDLGDV, from the C Terminus of the protein sequence according to **Immunogen**

NP 061934.2.

Peptide

C-LAKREDEFVDLGDV Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Dog

Reactive **Species**

Mouse

Human

Gene ID

54545

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png Grade

Western

Western blot: Approx 80-90kDa bands seen in Mouse Brain extracts. Recommended for use at 0.5-2µg/ml. **Blot**