

GOAT ANTI-MTMR12 ANTIBODY

SKU: EB05237



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	phosphatidylinositol-3-phosphate associated protein myotubularin related protein 12 3-phosphatase adapter subunit KIAA1682 protein KIAA1682 FLJ20476 PIP3AP 3-PAP phosphatidylinositol-3 phosphate 3-phosphatase adaptor subunit 3PAP MTMR12 protein
Accession ID	NP_061934.2
Immunogen	Peptide with sequence C-LAKREDEFVDLGDV, from the C Terminus of the protein sequence according to NP_061934.2.
Peptide Sequence	C-LAKREDEFVDLGDV
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Dog
Reactive Species	Mouse
Human Gene ID	54545
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 80-90kDa bands seen in Mouse Brain extracts. Recommended for use at 0.5-2µg/ml.