

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-MTMR3 ANTIBODY**

**SKU:** EB05193

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias

FYVE (Fab1 YGLO23 Vsp27 EEA1 domain) dual-specificity protein phosphatase|FYVE-DSP1|KIAA0371|myotubularin

related protein 3|MTMR3

**Names** 

Accession NP 066576

**Blocking** 

EBP05193 **Peptide** 

Immunogen Peptide with sequence C-DLELDKPIAATSN, from the C Terminus of the protein sequence according to NP 066576.

**Peptide** 

C-DLELDKPIAATSN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** Instructions

Refrigerated

**Predicted** 

Human, Mouse **Species** 

Human

8897 **Gene ID** 

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

Western Western blot: No signal obtained yet but low background observed in Mouse Brain and HeLa lysates at upto

Blot  $3\mu g/ml$ .