

GOAT ANTI-MTMR3 ANTIBODY

SKU: EB05193



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FYVE (Fab1 YGLO23 Vsp27 EEA1 domain) dual-specificity protein phosphatase FYVE-DSP1 KIAA0371 myotubularin related protein 3 MTMR3
Accession ID	NP_066576
Blocking Peptide	EBP05193
Immunogen	Peptide with sequence C-DLELDKPIAATSN, from the C Terminus of the protein sequence according to NP_066576.
Peptide Sequence	C-DLELDKPIAATSN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	8897
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Mouse Brain and HeLa lysates at upto 3µg/ml.